Urgotul®

INTERFACE LIPIDO-COLLOÏDE

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LIPIDO-COLLOID CONTACT LAYER

PRESENTATION

Urgotul is the product of a new and exclusive technology developed by URGO and called Technology Lipido-Colloid.

Urgotul is a non-occlusive dressing comprising of a polyester mesh impregnated with hydrocolloid particles dispersed in a petroleum jelly matrix. This mass contains carboxymethyl cellulose (hydrocolloid), paraffin, petroleum jelly and carrier polymers. **Urgotul** is presented in a sterile peelable sachet.

PROPERTIES

Urgotul is indicated for the management of acute wounds (burns, traumatic wounds, abrasions, post-operative wounds) and chronic wounds (leg ulcers, pressure ulcers and diabetic foot ulcers) at the granulation and epithelialisation stage. **Urgotul** is also intended for the management of epidermolysis bullosa. In contact with wound exudate, the hydrocolloid particles in the **Urgotul** mass interact with the petroleum jelly component to form a lipido-colloid interface.

This interface creates favourable, moist conditions for the healing process.

Due to a chemical composition which offers non-slip properties, but not greasy to touch, **Urgotul** does not adhere to either the wound or the surrounding skin. Removal of the dressing is painless and atraumatic for the wound.

METHOD OF USE

- Prepare the wound according to local protocol.
- Remove the protective transparent film.
- Apply **Urgotul** to the wound in a single layer.
- **Urgotul** can be cut with sterile equipment to fit the size of the wound.
- Cover **Urgotul** with a secondary dressing suitable to absorb the level of wound exudate.
- Hold secondary dressing in place with a conforming bandage, an adhesive tape or an elasticated tubular bandage.
- **Urgotul** can be left in-situ for up to 7 days depending upon the wound condition.
- For patients with epidermolysis bullosa, the dressing should be changed between 1 to 3 days.

PRECAUTIONS

- **Urgotul** adheres to latex surgical gloves. These should be moistened with normal saline to facilitate handling of **Urgotul**.
- If there are clinical signs of local infection, **Urgotul** may be used under the supervision of a health-care professional, in combination with suitable systemic antibiotics.

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